# ResistancePlus™ MG (550)


## For Applied Biosystems® (ABI) 7500 Fast

2000201

SpeedX Pty Ltd, Suite G16 National Innovation Centre  
 Australian Technology Park, 4 Cornwallis Street, Eveleigh,  
 Sydney, NSW 2015, Australia  
 Tel: +61 2 9209 4169, Email: [tech@speedx.com.au](mailto:tech@speedx.com.au)

MT Promedt Consulting GmbH, Altenhofstrasse 80,  
 66386 St. Ingbert, Germany  
 Tel: +49 6894 581020, Email: [info@mt-procons.com](mailto:info@mt-procons.com)

**FOR PROFESSIONAL USE ONLY\***

\*Not for sale in the USA

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### 1. Product description

The **ResistancePlus™** MG (550) kit is a qualitative multiplexed *in vitro* diagnostic real-time PCR test for the identification of *Mycoplasma genitalium* and detection of 5 mutations in the 23S rRNA gene (A2058G, A2059G, A2058T, A2058C, and A2059C, *Escherichia coli* numbering) that are associated with resistance to azithromycin (macrolide antibiotic). The **ResistancePlus™** MG (550) kit can be used on urine and urogenital swabs from male and female patients. It is intended to aid in the diagnosis of *M. genitalium* and detects mutations associated with azithromycin resistance in *M. genitalium*, and should be used in conjunction with clinical and other laboratory information.

Negative results do not preclude *M. genitalium* infections and do not provide confirmation of azithromycin susceptibility as there may be other mechanisms of treatment failure.

The **ResistancePlus™** MG (550) kit is intended to be used in professional settings such as hospitals, or reference or state laboratories. It is not intended for self-testing, home use, or point of care use.

### 2. Kit contents

Sufficient for 100 tests (20 µl)

Table 1. Kit contents		
Cap Colour	Content	Quantity
Brown	MG+23S mix, 20x	1 x 100 µl
White	Control Mix 2 <sup>^</sup> , 20x	1 x 100 µl
Orange	<b>Plex</b> Mastermix, 2x	1 x 1 ml
Blue	MgCl <sub>2</sub> , 50 mM	1 x 1.2 ml
Red	Internal Control Cells <sup>#</sup>	1 x 500 µl

<sup>^</sup> The Control Mix 2 contains a dye specific for the Applied Biosystems® (ABI) 7500 Fast  
<sup>#</sup> Store template tubes separately from oligo mixes, i.e. template or nucleic acid handling room

### 3. Instructions

Please refer to the **ResistancePlus™** MG kit Instructions for Use (IFU) for the detailed protocol (Ref No IF-IV0003). This document can be obtained by contacting [tech@speedx.com.au](mailto:tech@speedx.com.au).

It is recommended to use the **ResistancePlus™** MG Positive Control (Cat No 95001.) kit as positive controls for nucleic acid amplification, representing *M. genitalium* wild type and mutant *M. genitalium* with a 23S rRNA mutation (A2058G, A2059G, A2058T, A2058C, or A2059C *Escherichia coli* numbering). Please refer to the IFU for the detailed protocol.

### 4. Shipping and storage

The components of the **ResistancePlus™** MG (550) kit are shipped on dry ice or ice gel packs. All components should be stored at -20°C upon receipt. Excessive freeze/thawing is not recommended.

When stored under the recommended conditions and handled correctly, activity of the kit is retained until the expiry date stated on the label. Do not use past expiry date.

### 5. Interpretation of results

Please refer to the **ResistancePlus™** MG IFU for detailed interpretation of results.

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