
ResistancePlus[®] MG

For LightCycler[®] 480 Instrument II

20001L-01

MT Promedt Consulting GmbH, Altenhofstrasse 80,
66386 St. Ingbert, Germany
Tel: +49 6894 581020, Email: info@mt-procons.com

SpeedX Pty Ltd, Suite G16 National Innovation Centre
Australian Technology Park, 4 Cornwallis Street, Eveleigh,
Sydney, NSW 2015, Australia
Tel: +61 2 9209 4169, Email: tech@speedx.com.au

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1. Product description

The **ResistancePlus[®] MG** kit is a qualitative multiplexed in vitro diagnostic real-time PCR test for the identification of *Mycoplasma genitalium* and detection of 5 mutations in the 23S rRNA gene (A2058G, A2059G, A2058T, A2058C, and A2059C, *Escherichia coli* numbering) that are associated with resistance to azithromycin (macrolide antibiotic). The **ResistancePlus[®] MG** kit can be used on urine and urogenital swabs from male and female patients. It is intended to aid in the diagnosis of *M. genitalium* and detects mutations associated with azithromycin resistance in *M. genitalium*, and should be used in conjunction with clinical and other laboratory information.

Negative results do not preclude *M. genitalium* infections and do not provide confirmation of azithromycin susceptibility as there may be other mechanisms of treatment failure.

The **ResistancePlus[®] MG** kit is intended to be used in professional settings such as hospitals, or reference or state laboratories. It is not intended for self-testing, home use, or point of care use.

2. Kit contents

Sufficient for 100 tests (20 µl)

Table 1. Kit contents

Cap Colour	Content	Quantity
Brown	MG+23S mix, 20x	1 x 100 µl
White	Control Mix 1 [^] , 20x	1 x 100 µl
Orange	Plex Mastermix, 2x	1 x 1 ml
Blue	MgCl ₂ , 50 mM	1 x 1.2 ml
Red	Internal Control Cells [#]	1 x 500 µl

[^] The Control Mix 1 contains a dye specific for the LC480 II

[#] Store template tubes separately from oligo mixes, i.e. template or nucleic acid handling room

3. Instructions

Please refer to the **ResistancePlus[®] MG** kit Instructions for Use (IFU) for the detailed protocol (Ref No IF-IV0003). This document can be obtained by contacting tech@speedx.com.au.

The **PlexPCR[®] Colour Compensation** kit (Cat No 90001) **must** be run and applied before analysis of **PlexPCR[®]** kits on the LC480 II. This kit can be supplied on request. Please contact sales@speedx.com.au for more information.

It is recommended to use the **ResistancePlus[®] MG Positive Control** (Cat No 95001) kit as positive controls for nucleic acid amplification, representing *M. genitalium* wild type and mutant *M. genitalium* with a 23S rRNA mutation (A2058G, A2059G, A2058T, A2058C, or A2059C *Escherichia coli* numbering). Please refer to the IFU for the detailed protocol.

4. Shipping and storage

The components of the **ResistancePlus[®] MG** kit are shipped on dry ice or ice gel packs. All components should be stored at -20°C upon receipt. Excessive freeze/thawing is not recommended.

When stored under the recommended conditions and handled correctly, activity of the kit is retained until the expiry date stated on the label. Do not use past expiry date.

5. Interpretation of results

Please refer to the **ResistancePlus[®] MG** IFU for detailed interpretation of results.

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