
ResistancePlus[®] MG

ResistancePlus[®] MG for LightCycler[®] 480 Instrument II

2000125

MT Promedt Consulting GmbH, Altenhofstrasse 80,
66386 St. Ingbert, Germany
Tel: +49 6894 581020, Email: info@mt-procons.com

SpeedX Pty Ltd, Suite G16 National Innovation Centre
Australian Technology Park, 4 Cornwallis Street, Eveleigh,
Sydney, NSW 2015, Australia
Tel: +61 2 9209 4169, Email: tech@speedx.com.au

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Not for sale in the USA

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1. Intended Use

The **ResistancePlus[®] MG** kit is a qualitative real-time PCR test for the identification of *Mycoplasma genitalium* and detection of 5 mutations in the 23S rRNA gene (A2058G, A2059G, A2058T, A2058C, and A2059C, *Escherichia coli* numbering) that are associated with resistance to azithromycin (macrolide antibiotic). The **ResistancePlus[®] MG** kit is compatible with samples extracted on the instruments specified in the Instructions for Use (IFU) and real-time detection on the LightCycler[®] 480 Instrument II. The **ResistancePlus[®] MG** kit can be used with the specimen types listed in the Instructions for Use.

Negative results do not preclude *M. genitalium* infections and do not provide confirmation of azithromycin susceptibility as there may be other mechanisms of treatment failure.

The **ResistancePlus[®] MG** kit is intended to be used in professional settings such as hospitals, or reference or state laboratories. **It is not intended for self-testing, home use, or point of care use.**

2. Kit contents

Sufficient for 25 tests (20 µl)

Cap Colour	Contents	Quantity
Brown	MG+23S mix, 20x	1 x 25 µl
White	Control Mix 1, 20x	1 x 25 µl
Blue	Plex Mastermix, 2x	1 x 250 µl
Red	Internal Control Cells [#]	1 x 100 µl
Neutral	Nuclease Free Water	1 x 1 ml

[#] Store template tubes separately from oligo mixes, i.e. template or nucleic acid handling room

3. Instructions

Please refer to the **ResistancePlus[®] MG** kit Instructions for Use for the detailed protocol (Ref No IF-IV0003). This document can be obtained by contacting tech@speedx.com.au.

The **PlexPCR[®] Colour Compensation** kit (Cat No 90001) **must** be run and applied before analysis of **PlexPCR[®]** kits on the LC480 II. This kit can be supplied on request. Please contact sales@speedx.com.au for more information.

It is recommended to use the **ResistancePlus[®] MG Positive Control** (Cat No 95001) kit as positive controls for nucleic acid amplification, representing *M. genitalium* wild type and mutant *M. genitalium* with a 23S rRNA mutation (A2058G, A2059G, A2058T, A2058C, or A2059C *Escherichia coli* numbering). Please refer to the IFU for the detailed protocol.

4. Shipping and storage

The components of the **ResistancePlus[®] MG** kit are shipped on dry ice or ice gel packs. All components should be stored at -20°C upon receipt. It is recommended that freeze/thaw cycles are limited to 15.

When stored under the recommended conditions and handled correctly, activity of the kit is retained until the expiry date stated on the label. Do not use past expiry date.

5. Interpretation of results

Please refer to the **ResistancePlus[®] MG** IFU for detailed interpretation of results.

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