
PlexPCR[®] HSV-1&2, VZV

For LightCycler[®] 480 Instrument II

11101L-01

MT Promedt Consulting GmbH, Altenhofstrasse 80,
66386 St. Ingbert, Germany
Tel: +49 6894 581020, Email: info@mt-procons.com

SpeedX Pty Ltd, Suite G16 National Innovation Centre
Australian Technology Park, 4 Cornwallis Street, Eveleigh,
Sydney, NSW 2015, Australia
Tel: +61 2 9209 4169, Email: tech@speedx.com.au

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*Not for sale in the USA

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1. Product description

The **PlexPCR[®]** HSV-1&2, VZV kit simultaneously detects and differentiates Herpes simplex virus 1 (HSV-1), Herpes simplex virus 2 (HSV-2) and Varicella zoster virus (VZV) nucleic acids isolated and purified from cutaneous or mucocutaneous lesion swab specimens. It is a single well, multiplex real-time PCR (qPCR) assay and also contains a heterologous extrinsic control to monitor extraction efficiency and qPCR inhibition. The panel utilises **PlexZymes[®]** for specificity and superior multiplexing capability. The assay has been validated for use with samples extracted using the MagNA Pure 96 System (Roche) and real-time detection on the LightCycler[®] 480 Instrument II (LC480 II) (Roche).

2. Kit contents

Sufficient for 100 tests (20 µl)

Table 1. Kit contents

Cap Colour	Content	Quantity
Brown	HV3 Mix, 20x	1 x 100 µl
White	Control Mix 640 [†] , 20x	1 x 100 µl
Blue	Plex Mastermix, 2x	1 x 1 ml
Red	Internal Control Cells [‡]	1 x 500 µl
Neutral	Nuclease Free Water	1 x 1 ml

[†] The Control Mix 640 contains a dye specific for the LC480 II
[‡] Store template tubes separately from oligo mixes

3. Instructions

Please refer to the **PlexPCR[®]** HSV-1&2, VZV kit Instructions for Use (IFU) for the detailed protocol (Ref No IF-IV0002). This document can be obtained by contacting tech@speedx.com.au.

The **PlexPCR[®]** Colour Compensation kit (Cat No 90001) **must** be run and applied before analysis of **PlexPCR[®]** kits on the LC480 II. This kit can be supplied on request. Please contact sales@speedx.com.au for more information.

It is recommended to use Amplirun[®] DNA Controls (Vircell) as positive controls for HSV-1, HSV-2, and VZV (not supplied, Cat No MBC023, MBC024 and MBC048 respectively). Please refer to the IFU for the detailed protocol.

4. Shipping and storage

The components of the **PlexPCR[®]** HSV-1&2, VZV kit are shipped on dry ice or ice gel packs. All components should be stored at -20°C upon receipt. Excessive freeze/thawing is not recommended.

When stored under the recommended conditions and handled correctly, activity of the kit is retained until the expiry date stated on the label. Do not use past expiry date.

5. Interpretation of results

Please refer to the IFU for detailed interpretation of results.

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