PlexPCR® RespiVirus Positive Control

95002

MT Promedt Consulting GmbH

Altenhofstrasse 80

66386 St. Ingbert, Germany

Tel: +49 6894 581020, Email: info@mt-procons.com

SpeedX Pty Ltd

Suite G16, National Innovation Centre

Australian Technology Park, 4 Cornwallis Street, Eveleigh

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1. Intended use

The **PlexPCR**® RespiVirus Positive Control kit contains amplification controls. The positive controls contain synthetic DNA for Influenza A (FluA), Influenza B (FluB), Rhinoviruses (RhV), and Respiratory Syncytial Viruses A and B (RSV), Human metapneumovirus (HMPV), Adenoviruses B and C (AdV) and Human parainfluenza viruses 1, 2, 3 and 4 (HPIV 1-4) targets.

2. Kit contents

Cap colour	Contents [#]	Amplification control targets	Quantity
Orange	RV PC1	Flu A, Flu B, RhV, RSV A, HMPV, AdV B, HPIV 1	1 x 100 µl
Blue	RV PC2	RSV B, AdV C, HPIV 2	1 x 100 µl
White	RV PC3	HPIV 3	1 x 100 µl
Brown	RV PC4	HPIV 4	1 x 100 µl

[#] Store contents separately from mastermix setup area

3. Shipping and storage

The components of the **PlexPCR**® RespiVirus Positive Control kit are shipped on dry ice or ice gel packs. All components should be stored at -20°C upon receipt. It is recommended that freeze/thaw cycles are limited to less than 10.

When stored under the recommended conditions and handled correctly, activity of the kit is retained until the expiry date stated on the label. Do not use past expiry date.

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