



Evaluation of ResistancePlus GC to Inform Ciprofloxacin Treatment of Gonorrhoea

Litty Tan


Director of R&D

littyt@speedx.com.au

Evolution of Antimicrobial Resistance (AMR) in *Neisseria gonorrhoeae* (GC)



Antibiotics used for treatment



Emergence of resistance

WHO Global Action Plan for AMR GC



- 🧬 Increase awareness on correct antibiotic use
- 🧬 Effective prevention, diagnosis and control
- 🧬 Monitor treatment failures
- 🧬 Strengthen AMR surveillance
- 🧬 Effective drug regulations and prescription policies


Global action plan
to control the spread and impact of
antimicrobial resistance in *Neisseria gonorrhoeae*



- 🧬 Newer molecular methods for monitoring and detecting AMR
- 🧬 Alternative effective treatment regimens

Molecular detection of genetic markers for resistance/susceptibility could allow 'older' drugs to be used for gonorrhoea treatment

GC Susceptibility to Ciprofloxacin



% Susceptible to Ciprofloxacin




Heffernan H et al. (2015) Antimicrobial resistance and molecular epidemiology of gonococci in NZ, 2014-15
Kirkcaldy RD et al. (2016) MMWR Surveillance Summaries. July 15, 2016 / 65(7);1-19
Lahra MM et al. (2017) Australian Gonococcal Surveillance Programme Annual Report, 2017

GC Mechanisms of Fluoroquinolone Resistance



- Fluoroquinolones were predominant GC treatment in 1990s
 - 2007 - No longer recommended
- Quinolones act by inhibition of DNA gyrase and topoisomerase IV
 - DNA gyrase – heterotetramer of GyrA and GyrB
 - Topoisomerase IV – heterotetramer of ParC and ParE
- Bacteria develop resistance through mutations in the QRDR (quinolone resistance determining region)
 - GC: GyrA (aa positions 91, 95) and ParC (87)




GyrA Genotype Predicts Cipro Resistance in GC



Genomic Epidemiology in US, 2000-2013 isolates


Grad et al. JID 2016;214:1579-87



GyrA Genotype	Predictive Value
S91F	PPV of Ciprofloxacin Resistance = 98%
S91 WT	NPV of Ciprofloxacin Susceptibility = 99%

Genomic Survey in Europe, 2013 isolates

Harris et al. Lancet Infect Dis, published online May 15, 2018




GyrA Genotype Predicts Cipro Resistance in GC




Genomic Epidemiology in US, 2000-2013 isolates

Grad et al. JID 2016;214:1579-87



Genomic Survey in Europe, 2013 isolates

Harris et al. Lancet Infect Dis, published online May 15, 2018



GyrA S91/S91F is highly predictive of ciprofloxacin susceptibility/resistance

GC Molecular Diagnostics With *gyrA* Genotyping SpeedX Can Guide Ciprofloxacin Treatment



Advantages

- Oral treatment preferred
- Antibiotic stewardship of ceftriaxone
- Utilize existing drug while new drugs are still in development



Pooled estimate of real-time PCR *gyrA* genotyping tests for predicting GC susceptibility to ciprofloxacin

Allan-Blitz et al. Sex Transm Dis. 2017 May;44(5):261-265.

- Sensitivity: 98.2% (95% CI, 96.5–99.1%), Specificity: 98.6% (95% CI, 97.0–99.3%)




UCLA Health System has implemented a *gyrA* molecular assay to identify patients for ciprofloxacin treatment

Allan-Blitz et al. Sex Transm Dis. 2018 Apr;45(4):e18.

- 100% Ciprofloxacin cure rate for *gyrA* WT (n=25; 7 urethral, 7 pharyngeal, 7 rectal, 4 genital)

PlexPCR[®] Technical Advantages



Specific amplification and detection

Universal probe for easy multiplexing

Maintains sensitivity

Tan et al PLoS ONE (2017) and Patent: WO 2013/123552

Ideal molecular diagnostic technology for multiplexed detection of resistance/susceptibility markers

ResistancePlus GC



GC detection and ciprofloxacin resistance/susceptibility information

Rapid qPCR format (<1.5 hours)

Direct from clinical specimens

Specimen Types

Urine, Swabs

	Channel	Target
1 Well	1	<i>N. gonorrhoeae</i> (opa)
	2	<i>N. gonorrhoeae</i> (porA)
	3	gyrA S91 wild type
	4	gyrA S91F mutation
	5	Internal Control

ResistancePlus GC is a molecular diagnostic test for ciprofloxacin resistance/susceptibility

ResistancePlus GC (beta) Analytical Performance

Analytical sensitivity

- GC gyrA S91 wildtype – 15 geq/reaction
- GC gyrA S91F mutant – 15 geq/reaction

Analytical specificity

- 100% specificity:
 - *Neisseria spp.* isolates
 - Other organisms found in genital/throat/rectal sites

ResistancePlus GC (beta) on GC Clinical Isolates

 Australian isolates (top 70 most common genotypes from 2012)

		WGS	
		S91F mut	S91 WT
SpeedX	S91F mut	28	0
	S91 WT	0	42
	Total	28	42
Sensitivity		100.0% (95% CI 87.7-100.0%)	
Specificity		100.0% (95% CI 91.4-100.0%)	

		Ciprofloxacin AST	
		Resistant (R)	Susceptible (S)
SpeedX	S91F mut	27	1*
	S91 WT	0	42
	Total	27	43
Sensitivity		100.0% (95% CI 87.7-100.0%)	
Specificity		97.7% (95% CI 87.7-99.9%)	

* LS – less susceptible (S91F mutation by WGS)

High concordance to gyrA genotype & ciprofloxacin R/S phenotype

ResistancePlus GC (beta) Clinical Performance



- 🌀 University of Queensland Centre for Clinical Research (Brisbane, Australia)
- 🌀 416 specimens from 2017

	Urine	Genital	Throat	Rectal/ Anal	Other/ unknown	Total
Male	27	11	122	28	2	190
Female	18	184	22	1	1	226
Total	45	195	144	29	3	416

ResistancePlus GC (beta) Clinical Performance



 University of Queensland Centre for Clinical Research (Brisbane, Australia)

GC detection		GC Clinical Results	
		+	-
SpeedDx	+	123	1
	-	4	288
	Total	127	289
Sensitivity		96.9% (95% CI 92.1-99.1%)	
Specificity		99.7% (95% CI 98.1-100.0%)	

gyrA detection		In-house gyrA qPCR	
		Mutant	Wild type
SpeedDx	Mutant	20	0
	Wild type	0	70
	Indeterminate	0	1
	Total	20	71
Sensitivity		100.0% (95% CI 83.2-100.0%)	
Specificity		98.6% (95% CI 92.4-100.0%)	

Excellent clinical sensitivity and specificity

ResistancePlus GC – Future work



- Expected date for CE-IVD (Sept 2018)
- GRAND2 – GC clinical study

Potential Implementation Pathways



Screening or postal service

Asymptomatic
Test CT/GC



Confirm GC
ResistancePlus GC



GC + gyrA wild type
Treat with cipro



GC + gyrA mutant
Treat with ceft+azo

Partner testing

Partner test
ResistancePlus GC






GC + gyrA wild type
Treat with cipro



GC + gyrA mutant
Treat with ceft+azo

Acknowledgements

UQCCR, Australia

-  David Whiley
-  Lebogang Mhango
-  Ella Trembizki

SpeedX, Australia

-  Samantha Walker
-  Madeline Windsor
-  Simon Erskine
-  Elisa Mokany
-  Colin Denver

Come and visit SpeedX Booth #24