PlexPCR® Colour Compensation

LightCycler® 480 Instrument I & II and cobas® z 480 analyzer

REF 90001

EN **FOR PROFESSIONAL USE ONLY**

Not for sale in the USA

1. Intended use

The **PlexPCR®** Colour Compensation (CC) kit is for the generation of a **PlexPCR®** CC file on the LightCycler® 480 Instrument I & II (LC480 I & LC480 II) and cobas® z 480 analyzer (z 480) (Roche). All **PlexPCR®** kits and **ResistancePlus®** kits use a specific combination of dyes, and the **PlexPCR®** CC file **MUST** be applied for correct analysis of results.

2. Kit contents

96-well plate: Sufficient for 2 colour compensation files for 20 µl reaction volume
384-well plate: Sufficient for 5 colour compensation files for 10 µl reaction volume

Cap colour	Contents	Quantity
Blue	Plex Mastermix, 2x	1 x 1 ml
Brown	488 [†] CC mix, 20x	1 x 15 µl
Green	510 [†] CC mix, 20x	1 x 15 µl
Yellow	580 [†] CC mix, 20x	1 x 15 µl
Purple	610 [†] CC mix, 20x	1 x 15 µl
White	640 [†] CC mix, 20x	1 x 15 µl
Black	660 [†] CC mix, 20x	1 x 15 µl
Neutral	CC template [†]	1 x 400 µl
Light blue	Nuclease Free Water	1 x 1 ml

[†] Store template tubes separately from oligo mixes, i.e. template or nucleic acid handling room

[†] Emission channel on LC480 II (see Instructions for Use (IFU) for emission channel on LC480 I and z 480)

3. Instructions

Please refer to the **PlexPCR®** CC kit IFU for the detailed protocol (Ref No IF-IV0001). This document can be obtained by contacting tech@speedx.com.au.

The **PlexPCR®** CC file can be subsequently applied to all **PlexPCR®** or **ResistancePlus®** runfiles that have been created on the same instrument. Please refer to the **PlexPCR®** and **ResistancePlus®** kit IFU for the details on how to apply the **PlexPCR®** CC file.

4. Shipping and storage

The **PlexPCR®** CC kit is shipped on dry ice or ice gel packs. All components should be stored at -25°C to -15°C upon receipt.

When stored under the recommended conditions and handled correctly, activity of the kit is retained until the expiry date stated on the label. Do not use past expiry date.

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