PlexPCR® RespiVirus Positive Control

REF

95002

EN

FOR PROFESSIONAL USE ONLY

Not for sale in the USA

1. Intended use

The control reagent(s) consist of plasmid or synthetic RNA carrying target sequence(s) intended as a general-purpose non-assay specific reagent, for use as an amplification control in any qPCR reaction. This reagent does not contain any pathogen(s) and we do not make any claims, such as an assigned quantity (value) or specific quality parameters. This reagent has no diagnostic value and is provided as an accessory reagent, which possess no critical performance characteristics

2. Description

The **PlexPCR®** RespiVirus Positive Control kit contains amplification controls. The positive controls contain synthetic RNA for Influenza A (FluA), Influenza B (FluB), Rhinoviruses (RhV), and Respiratory Syncytial Viruses A and B (RSV), Human metapneumovirus (HMPV), Adenoviruses B and C (AdV) and Human parainfluenza viruses 1, 2, 3 and 4 (HPIV 1-4) targets.

3. Kit contents

Cap colour	Contents [#]	Amplification control targets	Quantity
Orange	RV PC1	Flu A, Flu B, RhV, RSV A, HMPV, AdV B, HPIV 1	1 x 100 µl
Blue	RV PC2	RSV B, AdV C, HPIV 2	1 x 100 µl
White	RV PC3	HPIV 3	1 x 100 µl
Brown	RV PC4	HPIV 4	1 x 100 µl

[#] Store contents separately from mastermix setup area

4. Shipping and storage

The components of the **PlexPCR®** RespiVirus Positive Control kit are shipped on dry ice or ice gel packs. All components should be stored at -25°C to -15°C upon receipt. It is recommended that freeze/thaw cycles are limited to less than 10.

When stored under the recommended conditions and handled correctly, activity of the kit is retained until the expiry date stated on the label. Do not use past expiry date.

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