# Positive Control Flu/RSV/SARS

**REF** 95006

EN

## FOR PROFESSIONAL USE ONLY

### 1 Intended Use

The Positive Control Flu/RSV/SARS swab is formulated with inactivated native Influenza A virus (FluA), Influenza B virus (FluB), Respiratory Syncytial Virus (RSV) and inactivated whole genome Severe Acute Respiratory Syndrome Coronavirus 2 (SARS-CoV-2) virus cDNA. The Positive Control Flu/RSV/SARS kit can be used as a full process control with respiratory virus nucleic acid detection assays, for use from sample extraction and purification to amplification and detection.

### 2 Kit contents

Packaging	Contents	Quantity
Aluminum foil pouch	Positive Control Flu/RSV/SARS	2 x 1 Swabs

### 3 Instructions

When using the Positive Control Flu/RSV/SARS in a test run, follow the same procedure for unknown specimens collected on a swab.

1. Remove swab from foil pouch.
2. To elute, completely submerge the swab end (tip) in Copan UTM media.
3. Snap the shaft end of the swab at the designated breakpoint, fasten screw cap and incubate for 1 min at room temperature.
4. Vortex the vial with the tip submerged for 30 seconds.
5. Process the liquid media from the eluted swab positive Control Flu/RSV/SARS for the nucleic acid extraction step.
6. After extraction, proceed with the nucleic acid test (NAT) by using the eluted nucleic acid test volume specified in the assay instruction for use. The positive control should be treated as though it were a clinical sample.

Once the foil pouch is opened, the Positive Control Flu/RSV/SARS swab should be used immediately.

### 4 Limitations

Protocols and results interpretation provided by manufacturers of test kits must be followed closely as deviations may produce unreliable results.

Failure to adhere to shipping and storage conditions may result in incorrect results.

### 5 Shipping and storage

The Positive Control Flu/RSV/SARS kit is shipped between 2°C to 30°C. The kit should be stored between 2°C to 30°C upon receipt.

When stored under the recommended conditions and handled correctly, performance of the kit is guaranteed until the expiry date stated on the label. Do not use past the expiry date.

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