

PlexPlus®

New
Innovative
Technology

Next Generation PCR Multiplexing



- ▶ Maximize PCR output with **PlexPlus®** technology that detects twice the targets compared to conventional PCR.
- ▶ **PlexPlus®** technology is achieved through specialized temperature-dependent probes, with fluorescence acquired at two temperatures across multiple wavelengths in real time, enabling affordable detection with higher sample throughput.
- ▶ The inaugural product, RespiV **PlexPlus®**, simultaneously detects 14 respiratory viral targets in a single well, and clearly discriminates between rhinovirus and enterovirus — whilst utilizing existing standard qPCR instrumentation in laboratories



RespiV *PlexPlus*^{®*}

One well. Double Detection. More Targets.

Assay Reagent RespiV *PlexPlus*^{®*} (Research Use Only)

Catalog Number 7300201

Description RespiV *PlexPlus*[®] is a 1-well, multiplex qPCR assay that amplifies nucleic acids of the following targets when present: influenza A virus (including H1N1 and H3N2 strains), influenza B virus (Yamagata and Victoria lineages), respiratory syncytial virus (subtypes A and B), severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2), human adenovirus (species B and C)[†], human parainfluenza virus (types 1 to 4), human metapneumovirus, rhinovirus (types A, B and C) and enterovirus (types A, B, C and D). This reagent includes an extraction/amplification internal control. The reagents are compatible with the Applied Biosystems[®] QuantStudio[™] 5 and 7 real-time PCR systems.

The RespiV *PlexPlus*[®] kit is for research use only (RUO) and is not intended for use in diagnostic procedures.

Sample Extraction Samples should be extracted as total nucleic acid.

Product Design				
Well	Channel	52 °C Temperature 1	76 °C Temperature 2	qPCR cycling duration
1	1	Influenza B virus	SARS-CoV-2	1 hr 55 minutes
	2	Rhinovirus/enterovirus	Enterovirus	
	3	Influenza A virus	Human adenovirus (B/C) [†]	
	4 [‡]	Respiratory syncytial virus (A/B)	Human parainfluenza virus (1/2/3/4)	
	5	Human metapneumovirus		
	6 [‡]	Internal control		
Instrument Compatibility				
Extraction		Liquid Handling	Amplification	
Magnetic bead-based		Lab self-validation	QuantStudio 5 and QuantStudio 7	
Workflow Flexibility				
Total Reactions		Reaction Volume		
100		20 µL		
Shipping & Storage				
Shipping			Storage	
Reagents — dry ice/ice packs			Store protected from light between -25°C to -15°C	

[†]Cross-detection of human adenovirus species A.

[‡]Custom calibration required. Refer to the relevant QuantStudio platform system manual for calibration instructions. Contact tech@speedx.com.au for assistance.

Americas - SpeedX Inc.
 +1 512 200 6918
 sales.us@speedx.com.au
 1400 Barbara Jordan Blvd,
 Austin TX 78723 USA

Australia - SpeedX Pty Ltd.
 +61 (0)2 9209 4170
 sales@speedx.com.au
 Suite 102, 4 Cornwallis Street,
 Eveleigh NSW 2015 Australia

Europe - SpeedX Ltd.
 +44 (0)330 445 0036
 sales.uk@speedx.com.au
 Acre House 11/15 William Road,
 London NW1 3ER United Kingdom