



# DermatoPlex\*

## Key target detection

**Test Reagent**

Dermato**Plex** (Research Use Only)

**Catalog Number**

7117192

**Description**

Dermato**Plex** is a 1-well real-time PCR multiplex for detection of common fungi causing dermatophytosis: **Trichophyton mentagrophytes complex** (including *T. mentagrophytes*, *T. interdigitale*, *T. indotineae*) and **Trichophyton rubrum complex** (including *T. rubrum*, *T. soudanense*, *T. violaceum*); and **Candida species** (including *C. albicans*, *C. glabrata*, *C. orthopsilosis*, *C. tropicalis* and *C. parapsilosis*).

The Dermato**Plex** kit is for research use only (RUO) and is not intended for use in diagnostic procedures.

Product Design			
Well	Channel	Target	qPCR cycling duration
1	1	<i>Trichophyton rubrum</i> complex ( <i>T. rubrum</i> , <i>T. soudanense</i> , <i>T. violaceum</i> )	1hr 15 minutes
	2	<i>Trichophyton mentagrophytes</i> complex ( <i>T. mentagrophytes</i> , <i>T. interdigitale</i> , <i>T. indotineae</i> )	
	3	Candida species ( <i>C. albicans</i> , <i>C. glabrata</i> , <i>C. orthopsilosis</i> , <i>C. tropicalis</i> and <i>C. parapsilosis</i> )	
	4	Internal Control	
Instrument Compatibility			
Extraction	Liquid Handling	Amplification	
Magnetic bead-based	Lab self-validation	Roche LC480 II	
		BioRad CFX96™ Dx and CFX96 Touch™	
		BioRad CFX384	
		ABI 7500 / 7500 Fast / 7500 Fast Dx	
Workflow Flexibility			
Total Reactions		Reaction Volume	
192		10 µL	
Shipping & Storage			
Shipping		Storage	
Reagents — dry ice/ice packs		Store protected from light between -25°C to -15°C	

**Australia - SpeedX Pty Ltd.**

+61 (0)2 9209 4170  
sales@speedx.com.au

Suite 102, 4 Cornwallis Street,  
Eveleigh NSW 2015 Australia

**Europe - SpeedX Ltd.**

+44 (0)330 445 0036  
sales.uk@speedx.com.au

Acre House 11/15 William Road,  
London NW1 3ER United Kingdom

**Americas - SpeedX Inc.**

+1 512 200 6918  
sales.us@speedx.com.au

1400 Barbara Jordan Blvd,  
Austin TX 78723 USA

\*Research Use Only  
Not for use in diagnostic procedures