

PlexPrime® GC+gyrA*

Test Reagent

PlexPrime® GC+gyrA (Research Use Only)

Description

Single-well multiplex for simultaneous detection of gene targets opa and porA (*N. gonorrhoeae* specific genes); gyrA S91F mutation; gyrA S91 wild type, and an internal control, using **PlexPCR®** technology

Product Design			
Well	Channel	Target	qPCR cycling duration
1	1	opa	1hr 15 minutes
	2	porA	
	3	gyrA S91F mutant sequence	
	4	gyrA S91 wild type	
	5	Internal Control	
Instrument Compatibility			
Extraction	Liquid Handling	Amplification	Cat No:
Magnetic bead-based	Lab self-validation	Roche LC480 II	Inquire
		CFX96 IVD / CFX96 Touch	Inquire
		ABI 7500 / 7500 Fast / 7500 Fast Dx / QuantStudio	7201031
Workflow Flexibility			
Total Reactions		Reaction Volume	
100		20 µL	
Shipping & Storage			
Shipping		Storage	
Reagents — dry ice/ice packs		Store protected from light between -25°C to -15°C	
		Instrument	Extraction
Sample Extraction		EZ1	EZ1 DSP Virus Kit
		MagNA Pure 96	MagNA Pure 96 DNA
		QIASymphony SP	DSP Virus/Pathogen Minikit

Australia - SpeedX Pty Ltd.

+61 (0)2 9209 4170
sales@speedx.com.au

Suite 102, 4 Cornwallis Street,
Eveleigh NSW 2015 Australia

Europe - SpeedX Ltd.

+44 (0)330 445 0036
sales.uk@speedx.com.au

Acre House 11/15 William Road
London NW1 3ER United Kingdom

Americas - SpeedX Inc.

+1 512 200 6918
sales.us@speedx.com.au

1400 Barbara Jordan Blvd
Austin, TX 78723, USA