


EntericPlex Parasite*

Sharpen Your Focus



Assay Reagent EntericPlex Parasite (Research Use Only)

Catalog Number 7151201

Description EntericPlex Parasite is a 1-well, multiplex qPCR assay that simultaneously detects and differentiates: *Entamoeba histolytica*, *Giardia intestinalis*, *Dientamoeba fragilis*, *Cryptosporidium* spp. and *Blastocystis* spp. This mix includes an internal control for extraction and amplification.

The EntericPlex Parasite kit is for research use only (RUO) and is not intended for use in diagnostic procedures.

Note: ENT-PAR Mix, 20x provided contains a fluorophore specific for the QuantStudio platforms.

Product Design			
Well	Channel	Target	qPCR cycling duration
1	1	<i>Entamoeba histolytica</i>	1 hr 40 minutes
	2	<i>Cryptosporidium</i> spp.	
	3	<i>Giardia intestinalis</i>	
	4†	<i>Dientamoeba fragilis</i>	
	5	<i>Blastocystis</i> spp.	
	6†	Internal Control	
Instrument Compatibility			
Extraction	Liquid Handling	Amplification	
Magnetic bead-based	Lab self-validation	QuantStudio 5 and QuantStudio 7	
Workflow Flexibility			
Total Reactions		Reaction Volume	
100		20 µL	
Shipping & Storage			
Shipping		Storage	
Reagents — dry ice/ice packs		Store protected from light between -25°C to -15°C	

†Custom calibration required. Refer to the relevant QuantStudio platform system manual for calibration instructions. Contact tech@speedx.com.au for assistance.

Americas - SpeedX Inc.
 +1 512 200 6918
 sales.us@speedx.com.au
 1400 Barbara Jordan Blvd,
 Austin TX 78723 USA

Australia - SpeedX Pty Ltd.
 +61 (0)2 9209 4170
 sales@speedx.com.au
 Suite 102, 4 Cornwallis Street,
 Eveleigh NSW 2015 Australia

Europe - SpeedX Ltd.
 +44 (0)330 445 0036
 sales.uk@speedx.com.au
 Acre House 11/15 William Road,
 London NW1 3ER United Kingdom